



PROSOMNUS [MOG]

The first precision Morning Occlusal Guide.

EVO™ [MOG]
Now
Available!

The ProSomnus [MOG] makes it easy to maintain or re-establish the bite in accordance with AADSM Guidelines⁽¹⁾ and select insurance reimbursement codes. Precision milled from digital records, the ProSomnus [MOG] is the only Morning Occlusal Guide that offers peace of mind that the bite can be recaptured for the duration of therapy – even if the initial Morning Occlusal Guide needs to be replaced.


Featuring two options for use with your sleep device cases based on clinical preference:

[MOG]: Design includes an anterior bite ramp.

[MOG] MIP: Design enables true, full occlusal contact.



[MOG]



[MOG] MIP

The objective of any [MOG], is to help re-establish the patients MIP after wearing a Mandibular Advancement Device. That can be done many ways. It usually involves a ramp that wedges the mandible back (just as their own teeth will) or it may have some soft material that helps the musculature seat the mandible back into place.

The precision, durability and reproducibility of the mandibular path and position are an advantage that provides peace of mind. You always have a guide to help re-establish the pre-treatment position.

Join the growing number of clinicians who trust ProSomnus devices for excellent patient experiences and outcomes.

CALL 844 537 5337 FOR A FREE STARTER KIT.



844 537 5337
ProSomnus.com
Patient Preferred OSA Therapy™



PROSOMNUS®
SLEEP TECHNOLOGIES