







# PROSOMNUS EVO™

For preferred OSA therapy.

Introducing ProSomnus EVO: The first OAT device to incorporate advanced materials, manufacturing robotics and artificial intelligence to advance the treatment of OSA.

- **100% of patients<sup>(1)</sup>** preferred EVO over previous devices and CPAP.
- **100% of providers<sup>(1)</sup>** would recommend EVO.
- **100% of providers<sup>(1)</sup>** stated EVO was easy to deliver.
- Preferred by patients, providers, physicians and payers.



**New MG6™ Technology** | ISO 10993 compliant | US Military drop test compliant | BPA Free  
Less porous, biocompatible and hygienic | Proprietary engineering and manufacturing

<sup>(1)</sup>ProSomnus Institutional Review Board Study data on file.

Join the growing number of clinicians who trust ProSomnus devices for excellent patient experiences and outcomes.







## CALL 844 537 5337 FOR A FREE STARTER KIT.



844 537 5337  
ProSomnus.com  
Patient Preferred OSA Therapy™



## Advancement Guide

5 DEVICES U0, U2,U4, L0, L1	COMBINATION	ADVANCEMENTS
	Upper 0 + Lower 0	0mm
	Upper 0 + Lower 1	1.0mm
	Upper 2 + Lower 0	2.0mm
	Upper 2 + Lower 1	3.0mm
	Upper 4 + Lower 0	4.0mm
	Upper 4 + Lower 1	5.0mm


When ordering additional advancement arches, please use the following guide for the next recommended arch sequence:

**U6** + L0 = 6.0mm

U6 + L1 = 7.0mm

**U8** + L0 = 8.0mm

U8 + L1 = 9.0mm



844 537 5337  
ProSomnus.com  
Patient Preferred OSA Therapy™



3 YEAR  
WARRANTY



**PROSOMNUS**  
SLEEP TECHNOLOGIES